
Comparison Charts Bacteria/Yeasts

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Ytor CFU/cm²

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
Comparison Chart Mould

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Diplides for industrial and Environmental Monitoring of Bacteria and Fungi

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Inoculation (Surfaces)

Remove cap, and withdraw combined cap/tongue unit. Do not touch culture gel. Holding the slide by the tabs at each end, press one side of the agar against the surface under test. The ribbed portion above the agar will bend to allow full contact.

Repeat with the second side on an area close to the first site.

Replace the slide in the tube, seal firmly, label as required, and incubate.

Inoculation (Fluids)

Remove cap/tongue unit from the tube, and without touching the gel, immerse the tongue in the fluid to be tested for about ten seconds, or hold the slide under running liquid or spray so that both sides are covered. Allow excess fluid to drain for a few seconds.

Replace in tube, seal firmly, label as required and incubate.

Incubation

Incubate upright as detailed for each product, **OR** at the temperature of the system if this is markedly different. Where other guidelines are laid down (eg. HSE) these should be followed.

Interpretation

Compare the incubated slide to the charts provided, reading the result from the closest matching picture. Fluids are calibrated in "Colony forming units" (CFU) per millilitre. Surface results are in CFU per square centimetre.

Bacteria and yeasts are both compared to the top chart. Fungi (Furry colonies) are compared to the lower chart.

Very high levels of infection may lead to a complete overgrowth on the slide. If this may be the case, comparing the slide to an unused example will show the difference.

Inkubering

Typ/kod Temp °C timmar

| Slide incubation | | | Agar 1 | | Agar 2 | |
|------------------|------|--------|----------|-----------------|-----------|--------------------|
| Code | Temp | Hrs | Colour | Growth | Colour | Growth |
| BT2 | 35C | 24-48 | Straw | Total Count | Straw | Total Count |
| BTM2 | 30C | 24-120 | Straw | Total Count | Brown | Yeasts & Moulds |
| PC2 | 35C | 24-48 | Straw | Total Count | Straw | Total Count |
| PCV | 35C | 24-48 | Straw | Total Count | Violet | Enterobacteriaceae |
| PO | 22C | 24-120 | Straw(*) | Total Count | Straw | Yeasts & Moulds |
| RBS | 30C | 24-120 | Straw | Total Count | Pink | Yeasts & Moulds |
| SC2 | 35C | 24-48 | Straw | Total Count | Red | Coliforms |
| SCN | 35C | 24-48 | Straw | Total Count | Red | Coliforms |
| SCT | 35C | 24-48 | Straw | Total Count | Red | Coliforms |
| TV | 35C | 24-48 | Straw | Total Count | Red/Brown | Enterobacteriaceae |
| BV | 35C | 24-48 | Yellow | Staphylococci | Violet | Enterobacteriaceae |
| BRB2 | 30C | 24-120 | Pink | Yeasts & Moulds | Pink | Yeasts & Moulds |
| TP | 35C | 24-48 | Straw(*) | Total Count | Straw | Pseudomonas |
| PC2T | 35C | 24-48 | Straw | Total Count | Straw | Total Count |
| PC2TN | 35C | 24-48 | Straw | Total Count | Straw | Total Count |
| BT3 | 35C | 24-48 | Straw | Total Count | Straw | Total Count |

Storage

Store in a cool dry place (Optimum 8-15 degrees c). Use before the expiry date on the box. Do not use slides if there is any sign of growth before opening. Slides remain usable while the agar is smooth and attached to the slide. A small amount of water may collect in the bottom of the tube. This is condensate from the agar and does not affect performance, provided the agar still fills the tray.

Disposal

Used slide should be incinerated, autoclaved or soaked in disinfectant for several hours before disposal.

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